User Manual

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
NON-CONTACT THERMOMETER

TERMÓMETRO SIN CONTACTO

English 02
Español 67

Contents

- 1 Thermometer

**ACTIV
ENERGY**

The contents include two AAA batteries.

Table of contents

Contents	2
General information	5
Reading and storing the user manual	5
Proper use	5
Contra-indications	6
Explanation of symbols.....	6
Safety	8
General safety information	8
Product description	18
Description of display	20
First use	22
Checking the thermometer and package contents	22
Information about body temperature	23
Start-up	25
Basic functions	27
Indicators for body temperature	27
Automatic deactivation.....	28
°Celsius/°Fahrenheit setting	28
Standby mode	29
Memory mode	29
Operation	31
Measuring forehead temperature.....	31
Measuring object temperature.....	33

Troubleshooting	36
Electromagnetic immunity (EMI).....	39
FCC Compliance Statement.....	39
Cleaning and maintenance	51
Cleaning	51
Calibration	54
Changing the battery	54
Storage	58
Technical data	59
Manufacturer information.....	64
Disposal.....	65
Disposing of the packaging.....	65
Disposing of the thermometer.....	65

General information

Reading and storing the user manual

This user manual accompanies this non-contact thermometer (referred to below only as the “thermometer”). It contains important information about start-up and handling.

Before using the thermometer, read the user manual carefully. This particularly applies to the safety instructions. Having the right measurement equipment is essential when it comes to obtaining accurate measurements. Failure to heed this user manual may result in severe injury or damage to the thermometer.

Store the user manual for future use. If you pass the thermometer on to third parties, please be absolutely sure to include this user manual.

Proper use

The thermometer is only designed for measuring the temperature of the human body as well as for measuring the temperature of liquids (e.g. milk or water) or surfaces. Body temperature measurements are taken on the forehead. The thermometer is intended for domestic and private use only, and not for commercial use.

The user is the proper operator of the thermometer. Only use the thermometer as described in this user manual. Any other use is deemed improper and may result in damage to property or even personal injury. The thermometer is not a children's toy.

The manufacturer or vendor accepts no liability for damage caused by improper or incorrect use.

Contra-indications

There are no known side effects linked to the use of the KFT509 thermometer.

Explanation of symbols

The following symbols and signal words are used in this user manual, on the thermometer or on the packaging.

This signal symbol/word designates a hazard with moderate degree of risk which may lead to death or severe injury if not avoided.

NOTICE!

This signal word warns against potential damages to property.

This symbol provides you with useful supplementary information about use.

Use and store in a dry state

Batch name:

LOT V4623KFT509

Serial number: **SN** 00001
(sequential number)

Type BF (body float) protection against electrical shock. A type BF applied part has greater protection against electric shock to the body, but is not suitable for direct cardiac application.

This symbol identifies information from the manufacturer (see chapter “Manufacturer information”).

2023-09-01

Production date:

2023-09-01 (year, month, day)

Safety

General safety information

Danger for children and persons with impaired physical, sensory or mental capacities (e.g. partially disabled persons, older persons with reduced physical and mental capacities) or lack of experience and knowledge (e.g. older children).

- Children or persons who lack knowledge and/or experience handling the thermometer or with impaired physical, sensory and/or mental abilities, must not use or play with the thermometer without supervision and instruction by a person who is responsible for their safety.
Attention! Risk of injury!

- Store the thermometer and accessories out of the reach of children and young persons under 18 years of age.
- Supervise children to prevent them from playing with the thermometer.
- Ensure that the thermometer does not end up in the hands of children. Children could ingest small parts or batteries and choke. Children could injure themselves when using the thermometer.
- Children's temperature measurements may only be taken by an adult.

WARNING!

Risk of injury!

Improper handling of the thermometer may result in injury.

- The thermometer must not be used in the event of a defect. Do not attempt to modify, disassemble or repair the thermometer or to replace parts of it.
- Check the thermometer for damage before each use. If there is an error or a defect (e.g. caused by children, pets or pests), the thermometer must not be used. In this case, contact the service center or the vendor.
- If you experience any irregularities (such as an allergic reaction) while using the thermometer, immediately stop using

- it and consult your physician.
- Do not engage in any other activity when performing a measurement.
 - If you are not sure how to use the thermometer, you should consult your physician before using it.
 - Keep this user manual for the entire service life of the product and include it if you pass the thermometer on to third parties. Make sure that the user manual is also accessible to third parties. The user manual is part of the thermometer. When passing it on to third parties, make sure they too are aware of the risks associated with the thermometer.
 - Avoid improper use and use that does not correspond to the designated purpose.

- Keep the infrared sensor clean and dry and avoid damages. This is the only way to obtain accurate measurements.
- This thermometer cannot serve as a substitute for medical advice or treatment! The measurements taken are only for reference. If you experience problems with your health, it is imperative that you consult your physician.
- Do not touch the infrared sensor and do not blow on it. Dirt on the infrared sensor could result in inaccurate measurement results.
- Do not use the thermometer if it is damaged. Using a damaged thermometer can result in injury, pose serious risks and lead to inaccurate measurement results.

- Wait a while before measuring temperature if the thermometer was stored in an area with a large temperature difference to the place where you intend to measure temperature. Both the thermometer and persons whose body temperature is being measured should have been in the same environment for at least 30 minutes.
- A greasy forehead can result in inaccurate measurements. You should therefore clean your forehead before taking your temperature. Then wait approx. 5–10 minutes before you perform a temperature measurement.
- Repeat the measurements approx. every 15 minutes if the readings are unusually low or if you feel they are inaccurate.

If possible, use a conventional clinical thermometer to take an oral or rectal temperature reading for comparison. Ask your physician if in doubt.

- To the extent possible, always measure the temperature of your body at the same point.

NOTICE!

Risk of damage!

Improper handling of the thermometer may result in damage.

- Do not use any accessories from other devices.
- Do not let the thermometer fall, do not shake it and do not expose it to shocks.
- Keep the infrared sensor of the thermometer clean and dry and avoid damages. This is the only way to obtain accurate measurements.

- Do not expose the thermometer to direct sunlight, fire, dirt, lint, dust, water, heat or extreme temperatures to prevent damage to the thermometer.
- Do not open the thermometer and do not attempt to repair it yourself in the event of disturbances or damage. This will void your warranty claim. The thermometer may only be repaired by authorized qualified professionals.
- Do not store the thermometer in areas where it is exposed to extremely high or low ambient temperatures (under $-4\text{ }^{\circ}\text{F}$ ($-20\text{ }^{\circ}\text{C}$) or over $131\text{ }^{\circ}\text{F}$ ($55\text{ }^{\circ}\text{C}$)) or in an extremely dry or moist environment as this can affect the accuracy of the measurements.

- Do not use the thermometer in a damp environment. Keep all liquids away from the thermometer during use and do not in any case immerse the thermometer in liquids.
- Do not use the thermometer if easily flammable substances and gases or explosives are nearby.
- During use, the thermometer may interfere with other electric or electrical equipment. Therefore, do not use it if other electric devices are nearby.
- Do not use a mobile phone near the thermometer when measuring temperature. Please note that portable and mobile HF (high-frequency) communication equipment (e.g. mobile phones) could interfere with medical electrical equipment.



- Do not use the thermometer at distances of less than 5 feet (1.5 meters) to a shortwave or microwave device or high-frequency surgical device.
- Do not use the thermometer in the mountains at altitudes of over 9,843 feet (3,000 meters).
- Medical electrical devices are subject to special precautionary measures with respect to the EMC (electromagnetic compatibility). Therefore, please observe the EMC notes concerning the installation and start-up of the thermometer.
- Remove the batteries from the thermometer if it is not being used for more than 3 months to prevent damage caused by leaky batteries.

Product description

The thermometer **1** is an infrared temperature-measuring device that you can use to easily perform precise measurements of body temperature on the forehead and of the temperature of liquids or surfaces. An acoustic signal indicates that the temperature measurement is complete; you can then read the temperature that has been measured on the LCD display **2**. The thermometer can save and display the data for the last 9 forehead temperature measurements.

The green or red illuminated display indicates whether the measured body temperature is in the normal range (green) or whether the body temperature is too high (red).


No special training or knowledge is required to use the thermometer.

- 2** LCD display
- 3** Button (⌘) for measuring temperature (including backlight)

- 4 Memory button (M)
- 5 Distance sensor
- 6 Infrared sensor

Description of display

- 7 Memory indicator (📄🔍)
(in memory mode)
- 8 Battery indicator (🔋)
- 9 Indicator for the room temperature (🏠)
- 10 Indicator for object temperature (🍼)

- 11 Segment indicator for the measured temperature
- 12 Indicator (☺) for measuring the temperature on the forehead

- In general, certain circumstances can affect body temperature such as age, clothing, outdoor temperature, physical activity, one's individual metabolism and time of day.
- Please note that certain conditions can affect measurements of body temperature taken on the forehead, for example, if your forehead was covered, after swimming or bathing or if your forehead has been subjected to particularly high or low temperatures.
- Holding the thermometer for a prolonged period of time may distort the measurement as it warms up the thermometer.

First use

Checking the thermometer and package contents

NOTICE!

Risk of damage!

If you are not cautious when opening the packaging with a sharp knife or other pointed objects, you may quickly damage the thermometer.

- Be very careful when opening it.
 1. Take the thermometer out of the packaging.
 2. Check to make sure that all parts are included (see **chapter “Contents”**).
 3. Check whether the thermometer or the individual parts exhibit damages. If this is the case, do not use the thermometer. Contact your branch store.

Information about body temperature

Body temperature refers to the internal temperature of the body (core temperature). It fluctuates over the course of the day and is lowest in the morning and highest in the late afternoon.

Normal body temperature is not reflected by a fixed value, but a range of values.

Age	Normal body temperature in °C	Normal body temperature in °F
0 to 2 years	36.4–38.0 °C	97.5–100.4 °F
3–10 years	36.1–37.8 °C	97.0–100 °F
11–65 years	35.9–37.6 °C	96.6–99.7 °F
over 65 years	35.8–37.5 °C	96.4–99.5 °F

The range for normal body temperature differs from person to persons and is affected by numerous factors: Age, gender, time of day, place and manner in which the temperature is measured, level of activity, mood and effect of medication and other factors. As a result, the body temperature of children is approx. 0.9 °F (0.5 °C) higher than that of adults.

Most people have a lower body temperature when they sleep; as the day progresses, their body temperature increases.

Measuring the temperature at different areas of the body will return different results: The body temperature measured rectally is normally approx. 0.54 °F (0.3 °C) higher than in the mouth and the temperature in the mouth is in turn 0.54 °F (0.3 °C) higher than the body temperature measured in the armpit.

To determine body temperature, it is recommended that you measure your body temperature when healthy on three to four consecutive days at different times of the day. This way, you can determine the normal body temperature, take note of it and, if necessary, compare the temperature readings with the normal readings taken for your body.


Start-up

- Before using the thermometer for the first time, carefully remove the protective film from the LCD display **2**.


The two batteries have already been inserted in the thermometer, but are protected by an insulating strip.

1. Before first use, open the battery compartment cover **13** on the back of the thermometer (see **Fig. 1**).


13 Battery compartment cover

2. Pull the insulating strip out and close the battery compartment cover again (see **Fig. 2**).

The thermometer will now emit 4 short signal tones to indicate that it has switched on and all indicators will initially appear on the LCD display for a short period (see **Fig. 3**).

The current room temperature will then appear on the display. The LCD display will also light up for approx. 15 seconds and the back-light for the button (8) for measuring temperature (3) will light up for approx. 10 seconds (see **Fig. 3**).

- If the thermometer is switched off, you can push any button. The back-light of the button for measuring temperature will then light up for approx. 15 seconds.
- If the thermometer is switched off, you can switch it on by pressing any button. The thermometer will then switch to standby mode. The room temperature will be shown on the LCD display in this mode.

- The thermometer will automatically switch to standby mode if there is no button activity within approx. 1 minute of taking the last temperature measurement. The thermometer will deactivate after another minute of inactivity.
- If “Lo” or “Hi” appears on the LCD display after taking a measurement, the temperature measured has exceeded or underrun the permissible temperature range.

Basic functions

Indicators for body temperature

The button (8) for measuring temperature has a backlight, which immediately illuminates in different colors after measuring the body temperature to indicate the range of your body temperature.

Green backlight:

The body temperature measured is below 100.3 °F (37.9 °C).

Red backlight:

The body temperature measured is within the range of an elevated temperature of 100.4 °F (38 °C) to 109.3 °F (42.9 °C).

If the thermometer is switched off, you can push any button. The backlight of the button for measuring temperature will thus be activated and will light up in blue for approx. 15 seconds.

Automatic deactivation

By default, the thermometer is set to switch off approx. 60–120 seconds after the last button activity.

°Celsius/°Fahrenheit setting

1. Push any button to switch to standby mode.
2. Hold the memory button (**M**)  down for approx. 2 seconds until the temperature measurement/indicator switches from °C (degrees Celsius) to °F (degrees Fahrenheit).


A short signal tone will sound to indicate you have switched from °C to °F.

3. Follow exactly the same steps to switch from °F (degrees Fahrenheit) to °C (degrees Celsius).

Two short signal tones will sound to indicate you have switched from °F to °C.

Standby mode

The room temperature indicator (🏠) **9** will appear on the LCD display **2** in standby mode (see **Fig. 1**).


Memory mode


The thermometer will save the last 9 body temperature measurements. The current measurement will appear as the last result in the memory. If all 9 memory slots are occupied, the oldest measurement will be deleted from the memory.

1. In standby mode, push the memory button (M) **4**.

The memory indicator (MEM) **7** and the number of the last measurement saved (e.g. 6) will now appear shortly on the LCD display **2** (see **Fig. 2**).

The measured temperature together with the indicator (☺) for measuring the temperature on the forehead **12** will then appear (see **Fig. 3**).

2. Each time the memory button (**M**) is pushed, first the number and then the data saved for the next memory slot will be shown shortly in descending order.

If you push and hold the memory button (**M**) for approx. 2 seconds while the saved measurement data is shown, the temperature in the memory will switch from °C (degrees Celsius) to °F (degrees Fahrenheit) or vice-versa. The switch from °C to °F will be indicated by a short signal tone; two signal tones will sound to indicate the switch from °F to °C.

Operation

There is no biological risk or hazard to the user if the thermometer comes into contact with the skin during use.


- Before use, let the thermometer adjust to the temperature of the room or surroundings to prevent incorrect measurement results.

Measuring forehead temperature


If possible, always take the temperature measurement in the same way for consistent measurement results.


1. Before taking a measurement, make sure that the infrared sensor **6** is clean.
2. Hold the thermometer so that your thumb is on the button **(8)** for measuring temperature **3** and push the button once to switch the thermometer on (see **Fig. 1**).

3. Hold the thermometer **1** in the middle of your forehead (see **Fig. 2**).


4. Hold the thermometer **1** steady and then push the button **(8)** for measuring temperature (see **Fig. 3**).

The distance sensor **5** will then measure the distance between the thermometer and your forehead. The ideal distance amounts to 0–1.97" (0–5 cm).

A continuous signal will sound if the distance is too great. If the distance is correct, a signal tone will sound twice and the temperature will be measured.

The indicator **(☺)** for measuring the temperature on the forehead **12** and the measured temperature will appear on the LCD display **2** (see **Fig. 4**).

At the same time, the backlight for the button (8) for measuring temperature will illuminate, which indicates the range for the body temperature measured (see chapter “Indicators for body temperature”).


If you want to immediately switch back to standby mode, you can briefly push the memory button (M). The thermometer will also switch back to standby mode if no button is pushed for approx. 1 minute.

Measuring object temperature


Please note: This measurement method is not suitable for measuring body temperature.

1. Before taking a measurement of an object (liquid/surface), make sure that the infrared sensor (6) is clean.
2. Hold the thermometer so that your thumb is on the button (8) for measuring temperature (3) and push the button once to switch the thermometer on (see **Fig. 1**).

3. Push the button (8) for measuring temperature and the memory button (M) (4) at the same time. The indicator for object temperature (10) will appear together with the indicator “Food” on the LCD display (2) (see **Fig. 2**).


4. Hold the thermometer (1) with the infrared sensor so that it can measure the temperature of the object within a distance of 0–1.97" (0–5 cm) (see **Fig. 3**).

5. Now hold down the button (8) for measuring temperature. A signal tone will indicate that the measurement has started.
-

6. Hold the thermometer at the same location for a few seconds until the measurement results have stabilized.

Once the temperature measurement is complete, the object temperature measured will appear on the LCD display **2** (see **Fig. 4**).


7. You can now let go of the button (8) for measuring temperature.

If you want to immediately switch back to standby mode, you can briefly push the memory button (**M**). The thermometer will also switch back to standby mode if no button is pushed for approx. 1 minute.

Troubleshooting

Fault	Cause	Solution
<p>The batteries have been inserted, but the LCD display  is blank.</p>	<p>There is debris in the battery compartment.</p>	<p>Remove the debris from the battery compartment.</p>
	<p>The batteries are empty or were inserted with incorrect polarity.</p>	<p>Replace the batteries with full batteries. Make sure the polarity is correct!</p>
	<p>The battery contacts are not aligned with the terminals.</p>	<p>Check the battery contacts.</p>
	<p>There is a fault in the electronics.</p>	<p>Remove the batteries and put them back in after approx. 3 seconds.</p>
	<p>The temperature measured is less than 89.6 °F/32.0 °C or higher than 109.3 °F/42.9 °C.</p>	<p>Observe the notices in the user manual.</p>

Fault	Cause	Solution
 	The thermometer is defective (hardware error).	Contact the vendor/distributor.
	The operating temperature is outside of the valid range.	Observe the operating conditions.
The temperature measured is not accurate or you doubt the accuracy of the measurement.	Check whether the infrared sensor 6 is dirty.	Clean the infrared sensor (see chapter “Cleaning and maintenance”).
	Check whether the thermometer is properly used.	Observe the notices on operating the thermometer and obtaining a correct measurement.

Troubleshooting

Fault	Cause	Solution
	Check if the thermometer and patient have been in the same room for at least 30 minutes before performing the measurement.	Make sure that the thermometer and patient have been in the same room for at least 30 minutes before performing the measurement.
	The thermometer is being used outdoors.	Only use the thermometer in indoor areas.
	You have held the thermometer in your hands for a prolonged period.	Place the thermometer e.g. a table so that it can adjust to the room temperature.
The thermometer behaves unexpectedly during operation.	The thermometer is defective.	If any irregularity occurs while the thermometer is being used, report this to the service center or vendor.

Electromagnetic immunity (EMI)

FCC Compliance Statement

Responsible party

Company name: Greater Concepts by
Design

Address: 27825 Ellie Oak Lane,
Spring Texas 77386, USA

Phone number: 877-944-3466

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to the product not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates,

uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.
- This device should not be used next to other devices; nor should it be placed on top of other devices. If contact cannot be avoided, this device should be monitored to ensure it operates properly in the arrangement used.
- Using accessories other than those specified or provided by this appliance's manufacturer may result in increased electromagnetic emissions or reduced

electromagnetic immunity of this device as well as in improper operation.

- Portable RF communications equipment (including peripheral components such as antenna cables and external antennas) should be placed at least 11.8 inches (30 cm) away from the device's parts. This also applies to the cable specified by the manufacturer. Otherwise, there may be a drop in the performance of this device.
- Strong electromagnetic interference usually occurs when there is a relatively dry operating environment. This can have the following effects on the device:
 - the device stops emitting;
 - the device turns off;
 - the device restarts;

The aforementioned phenomenon does not have any effect on the basic safety and performance of the device. The user can use it as described in the user manual. If you wish to prevent the aforementioned phenomenon, use the device in the environment specified in the user manual.

Table 1

Explanation – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.		
Emission test	Conformity	Electromagnetic environment – instructions
RF emissions CISPR 11	Group 1	The device only uses RF energy for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Electromagnetic immunity (EMI)

Emission test	Conformity	Electromagnetic environment – instructions
RF emissions CISPR 11	Class B	The device is suitable for use in residential areas and in areas directly connected to the public low-voltage power supply that supplies power to residential buildings.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ Flickering IEC 61000-3-3	Not applicable	

Table 2

Explanation – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – instructions
Electrostatic discharge IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If the floors are covered with synthetic material, the relative humidity should be at least 30 %.
Fast transient electrical disturbances/bursts IEC 61000-4-4	± 2 kV for power supply cables ± 1 kV for input/output cables	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.

Electromagnetic immunity (EMI)

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – instructions
Surge IEC 61000-4-5	<p>± 0.5 kV, ± 1 kV cable(s) to cables</p> <p>± 0.5 kV, ± 1 kV, ± 2 kV cable(s) to cables</p>	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<p>0 % UT; half cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°</p> <p>0 % UT; 1 cycle and 70 % UT; 25/30 cycles</p> <p>Single phase: at 0°</p> <p>0 % UT; 250/300 cycles</p>	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user must continue to use the device during power outages, it is recommended that the device be powered by an uninterruptible power supply or a battery.

Electromagnetic immunity (EMI)

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – instructions
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Important: UT is the (alternating) supply/line voltage prior to application of the test levels.			

Table 3

Explanation – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – instructions

Electromagnetic immunity (EMI)

<p>Conduct- ed RF IEC 61000-4-6</p>	<p>3 V 0.15 MHz to 80 MHz 6 V in ISM and ama- teur radio bands between 0.15 MHz and 80 MHz</p>	<p>Not ap- plicable</p>	<p>Portable and mobile RF com- munications equipment should be used no closer to any part of the device than the recom- mended separa- tion distance calculated from the equation applicable to the frequency of the transmitter.</p>
<p>Radiated RF IEC 61000-4-3</p>	<p>10 V/m 800 MHz to 2.7 GHz</p>	<p>10 V/m</p>	<p>Recommen- ded separation distance $d = 1.2 \sqrt{P}$ 150 kHz to 80 MHz $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum out- put power rating of the transmitter in watts (W) according to the transmitter man- ufacturer and d is the recommend- ed separation distance in me- ters (m).</p>

Electromagnetic immunity (EMI)

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – instructions
			<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Disturbances are possible within the vicinity of devices with the following symbol:</p> 
<p>NOTICE 1 In the case of 80 MHz and 800 MHz, the higher frequency range shall apply.</p> <p>NOTICE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

a. The field strength of stationary transmitters such as base stations for radio telephones (mobile/cordless phones) and land mobile radio services, amateur stations, AM and FM radio and television stations cannot be precisely predetermined in a theoretical respect. An electromagnetic survey of the site should be considered to evaluate the electromagnetic environment based on RF transmitters. If the field strength measured at the location where the device is used exceeds the aforementioned valid RF compliance level, the device should be monitored to ensure proper operation. If you observe unusual operations, additional measures may be necessary, for example, readjusting or moving the device.

b. Over the frequency range of 0.15 MHz to 80 MHz, the field strengths should amount to less than 3 V/m.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device

The device model is intended for use in an electromagnetic environment in which RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device model as recommended below, according to the maximum output power of the communications equipment.

Electromagnetic immunity (EMI)

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	$d = 1.2 \sqrt{P}$ 0.15 MHz to 80 MHz	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz	$d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTICE 1 In the case of 80 MHz and 800 MHz, the higher frequency range shall apply.

NOTICE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Cleaning and maintenance

Cleaning

NOTICE!

Risk of short circuit!

Water or other liquids that have penetrated the housing may cause a short circuit.

- Never submerge the thermometer in water or other liquids.
- Make sure that no water or other liquids penetrate the housing.

NOTICE!

Risk of damage!


Improper handling of the thermometer may result in damage to the thermometer.

- Do not use any aggressive cleaners, brushes with metal or nylon bristles or sharp or metallic cleaning utensils such as knives, metal sponges and the like. They could damage the surfaces.
- Never attempt to clean the thermometer with a scouring agent, gasoline or a solvent.


- For accurate measurements, it is very important that the infrared sensor of the thermometer is clean and free of scratches.
- Clean and disinfect the thermometer before handing it over to other persons, to prevent diseases from being transmitted.

1. Before taking every measurement, check the infrared sensor **6** and, if necessary, clean it with a cotton swab dipped in a small amount of medical alcohol and directly afterwards with a clean, dry cotton swab (for about 10 seconds with 75 % alcohol) (see **Fig. 1**).

2. Carefully clean the surfaces of the thermometer **1** (and of accessories) at regular intervals (e.g. after every 20th use) with a soft, damp cloth and a small amount of rubbing alcohol (see **Fig. 2**). Please ensure that no moisture enters the thermometer. A mild cleaning agent may be used for stubborn dirt. Clean the thermometer only when it is switched off. Remove the batteries before cleaning the thermometer.

3. Allow the thermometer to dry thoroughly. Do not use chemical cleaners or solvents to clean the thermometer.

Calibration

The thermometer was calibrated in the factory following manufacture. If the thermometer is used in accordance with the provisions and notices in the user manual for private and non-commercial use, there is no need to calibrate it at regular intervals. If, however, you doubt the accuracy of the measurements, you should immediately contact your vendor.

- Never attempt to disassemble the thermometer or make changes as this could result in risks and in such a case, also void any warranty claim.

Changing the battery

Danger of explosion!

Improper handling of the batteries may cause an explosion.

- Only replace the batteries with the same or an equivalent battery type.

- Do not expose the batteries to excessive heat such as from sunshine or fire.
- Protect the batteries from mechanical damage.
- The enclosed batteries must not be charged or reactivated with other means, they must not be dismantled, disposed of in fire or short circuited.

WARNING!

Risk of injury!

Handling batteries improperly may result in injury.

- Store the batteries so that they are not accessible for children.
- If swallowed, batteries could pose a risk of death. Therefore, store the batteries so that they are not accessible for small children. If a battery

is swallowed, you must seek immediate medical attention.

- If battery acid has leaked, avoid contact with your skin, eyes and mucous membranes. In the event of contact with acid, immediately flush the applicable areas with plenty of clean water and promptly consult a physician.

NOTICE!

Risk of damage!

Improper handling of the batteries may result in damage.

- Make sure the polarity is correct when inserting the batteries.
- Only insert batteries of the same type; do not use different types of batteries or old batteries with new batteries.


- Only use type AAA alkaline batteries. Do not use any rechargeable batteries.
- Remove the old batteries from the battery compartment if they are empty or if you do not intend to use the thermometer for a prolonged period of time.

If the battery symbol () appears on the LCD display **2** or the thermometer no longer functions, the batteries should be replaced immediately.

All settings and entries in the memory will be erased if the batteries are empty or when changing weak batteries.

1. To open the battery compartment cover **13**, push against the locking clip on the back of the thermometer **1** and take the battery compartment cover off (see **Fig. 1**).

2. Remove the empty batteries.
3. Insert two new alkaline batteries (type AAA). When inserting the batteries, make sure that the polarity is correct (see marking/stamped information in the battery compartment).
4. Put the battery compartment cover back in and push it against the thermometer housing until it locks back into place.

Storage

If you do not use the thermometer for a prolonged period of time, please observe the following notes:

- The thermometer requires no servicing.
- Take the batteries out of the battery compartment if you do not use the thermometer for more than 3 months.
- Do not store the thermometer near heat sources.
- Keep the thermometer away from children, pets and pests.
- Do not expose the thermometer to direct sunlight, fire, dirt, lint, dust, water, heat or extreme temperatures to prevent damage to the thermometer.
- Note the information regarding storage/transport conditions in the chapter “Technical data”.
- Do not place any heavy objects on the thermometer.

- Do not allow the thermometer to fall, do not expose it to shocks and do not throw it.

Technical data

Model/type:	KFT509
Product code:	829852
Dimensions (L × W × H):	approx. 5.12 × 1.77 × 1.98" (130.1 × 45 × 50.3 mm)
Weight:	approx. 0.14 lb (62 g) (without batteries)
Material:	Plastics, metals
Power supply:	DC 3V, 2× AAA batteries (V = volts)
IP-Code:	IP22 Meaning of the first digit: Protected against solid foreign bodies with a diameter greater than 0.5" (12.5 mm). Protected against access with a finger. Second digit: Protection against dripping water falling vertically when the housing is tilted up to 15°.

Technical data

Measurement method:	contactless infrared measurement
Application part:	Measuring sensor in the upper device area
Measurement time:	approx. 2 seconds
Measuring distance:	approx. 0–1.97" (0–5 cm)
Display:	LCD display, indicator shown in °C (degrees Celsius) or °F (degrees Fahrenheit), resolution: 0.1 °C/°F
Memory:	Storage of the last 9 body temperature measurements
Software version:	V1
Energy-saving function:	Depending on the enabled mode, the thermometer automatically deactivates after approx. 60 seconds.

Operating life of device: Approx. three years (with 10 measurements per day)

Operating life of the battery: approx. 6 months with 5 measurements per day

Reference body site: Core

Operating mode: Adjusted mode

Body temperature:

Measurement range: 89.6 °F–109.3 °F
(32.0 °C–42.9 °C)

Laboratory accuracy: from 95.0 °F–107.6 °F
(35.0 °C–42.0 °C):
+/-0.4 °F (+/-0.2 °C)
other ranges:
+/-0.5 °F (+/-0.3 °C)

Object temperature:

Measurement

range: -7.6 °F–176.0 °F
(-22.0 °C–80.0 °C)

Measurement

accuracy: +/-4.0 °F (+/-2.0 °C)

Room temperature:

Measurement

range: 41.0 °F–139.8 °F
(5.0 °C–59.9 °C)

Measurement

accuracy: +/-2.0 °F (+/-1.0 °C)

Operating conditions:

Atmospheric

pressure: 700 hPa–1,060 hPa

Forehead and tar-

get measurement: 59.0 °F–104.0 °F
(15.0 °C–40.0 °C) at
15 %–85 % relative
humidity

Clinical repeatability:


Measuring unit	Group A1	Group A2	Group B	Group C
Clinical distortion Δ_{cb}	-0.185	-0.175	-0.173	-0.182
Measurement uncertainty L_A	0.357	0.349	0.408	0.391
Clinical repeatability σ_r	0.102	0.081	0.118	0.212

Age group

Age

A1	0–3 months
A2	3 months–1 year
B	>1 year–<5 years
C	>5 years

Storage/transport conditions:

Storage/transport temperature:

-4 °F–131 °F
(-20 °C–55 °C) at
15 %–90 % relative
humidity

Atmospheric pressure:
700 hPa–1,060 hPa

This infrared thermometer meets the ear canal requirements established in ASTM Standard E 1965-98 for the thermometer system. Full responsibility for the conformance of the

Manufacturer information

product to the standard is assumed by Shenzhen Dongdixin Technology Co., Ltd., Floor 1-2, No.3 Building, Fanshen Xusheng Industrial Estate Xilixiaobaimang 518108 Nanshan District, Shenzhen P. R. China. ASTM laboratory accuracy requirements for the thermometer only in the display range of 98 °F to 102.2 °F (37 °C to 39 °C) for infrared thermometers is ± 0.4 °F (± 0.2 °C), whereas for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E 667-86 and E 1112-86 is ± 0.2 °F (± 0.1 °C).

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Disposal

Disposing of the packaging

- Dispose of the packaging separated into single type materials. Dispose of paper-board and cardboard with waste paper and plastics with recyclable waste.

Disposing of the thermometer

- Should the thermometer no longer be capable of being used at some point in time, **dispose of it in accordance with the regulations in force in your city or state.**
- Please ensure your recycling information applies to local regulations and the EPA recommendations (www.epa.gov).

Batteries and rechargeable batteries may not be disposed of with household waste!

As the end user you are required by law to bring all batteries and rechargeable batteries, regardless whether they contain harmful substances* or not, to a collection point run by the city, county, or state, or to a retailer, so that they can be disposed of in an environmentally friendly manner.

*labeled with: Cd = cadmium, Hg = mercury, Pb = lead

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Shenzhen Dongdixin Technology Co., Ltd., Floor 1-2, No.3 Building, Fanshen Xusheng Industrial Estate, Xilixiaobaimang, 518108 Nanshan District, Shenzhen, República Popular China

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US

Made in China

MANUFACTURER:
SHENZHEN DONGDIXIN TECHNOLOGY CO., LTD.
FLOOR 1-2, NO.3 BUILDING,
FANSHEN XUSHENG INDUSTRIAL ESTATE
XILIXIAOBAIMANG 518108 NANSHAN DISTRICT,
SHENZHEN
P. R. CHINA

DIST. & SOLD EXCLUSIVELY BY:
ALDI
BATAVIA, IL 60510

AFTER SALES SUPPORT
SERVICIO POSVENTA

829852

US

877-944-3466

greaterconcepts@olive.com

MODEL/MODELO:
KFT509

01/2024

3

YEARS WARRANTY
AÑOS DE GARANTÍA